



CPAP Resupply Schedule

CPAP supplies are made to be disposable and are manufactured with replacement schedules in mind, with proven medical reasoning behind the industry standard replacement guidelines. Cleaning supplies daily can help them last their longest and decrease the likelihood of bacteria forming, but adhering to a resupply schedule is vital to comfortable and successful CPAP therapy.

Medicare and many insurance plans will provide coverage for CPAP supplies after a deductible has been met, and Aeroflow Healthcare can help patients obtain their CPAP supplies on schedule.

CPAP Replacement Reference Guide:

Disposable Filters	2 Times Every Month
Nasal Cushions and Pillows	2 Times Every Month
Full Face Cushions	1 Time Every Month
Oral, Nasal, Nasal Pillow, Full Face Masks	1 Time Every 3 Months
Tubing	1 Time Every 3 Months
Headgear and Chinstraps	1 Time Every 6 Months
Non-Disposable Filters	1 Time Every 6 Months
Humidifier Chamber	1 Time Every 6 Months
CPAP Machine	1 Time Every 5 Years

If you have questions about cleaning or replacing your supplies, contact an Aeroflow Healthcare representative at 888-345-1780 or visit our cash-pay resupply website at CheapCPAPSupplies.com.